



- Keep clean when not using
- Clean and store with the bristles straight (with the exception of the stippling brush, which is meant to be frayed to replicate tissue texture)
- Clean with water at end of day and store as mentioned above
- Don't let materials dry to a solid state while in the brush bristles
- Do not clean with any solvent such as alcohol, mineral spirits, turpentine or acetone
- Use the Spatula to mix MiYO Colors and Glaze in the jar (as shown below), do NOT use the brushes to stir



Jensen Dental
miyoworld.com | jensendental.com
800.243.2000

