

How to identify an Ivoclar Programat CS model



Ivoclar Programat Furnace Model CS
Model noted above screen, top right of base

Crystallization	Vacuum	Sintering	Open Programs
Yes	Yes	No	Yes



Ivoclar Programat Furnace Model CS2
Model noted above screen, top left of base

Crystallization	Vacuum	Sintering	Open Programs
Yes	Yes	No	Yes



Ivoclar Programat Furnace Model CS3
Model noted above screen, top left of base

Crystallization	Vacuum	Sintering	Open Programs
Yes	Yes	No	Yes



Ivoclar Programat Furnace Model CS4
Model noted top center of hood

Crystallization	Vacuum*	Sintering	Open Programs
Yes	No	Yes	Yes

*no vacuum, MiYO Structure NOT recommended



Ivoclar Programat Furnace Model CS6
Model noted top center of hood

Crystallization	Vacuum	Sintering	Open Programs
Yes	Yes	Yes	Yes

*singles/small cases only due to muffle size and overall furnace design