


[NEW]



Shade conversion. Simplified.

Out of the jar shade conversion made easy.

-  No complicated formulas to mix
 -  Highly esthetic end results
 -  Simple application process
-

Convert monolithic VITA® classic to select 3D-MASTER® shades, straight from the jar.

Easy as 1...2...miYO^{3D}

11 Translucent Colors designed to achieve 12 VITA 3D-MASTER® Shades.



1

SELECT

Target VITA 3D-MASTER® Shade
1M1
1M2
2M1
2M2
2M3
3M1
3M2
3M3
4M1
4M2
4M3
5M3

2

CHOOSE

Starting Substructure Shade
BL3
B1
BL3
A1
A1
BL3
C1
A1
C1
C2
A3.5
A4

3

APPLY

MiYO® 3D Trans Shade (1-11)
MiYO 3D Trans Shade 1
MiYO 3D Trans Shade 2
MiYO 3D Trans Shade 3
MiYO 3D Trans Shade 4
MiYO 3D Trans Shade 5
MiYO 3D Trans Shade 6
MiYO 3D Trans Shade 2
MiYO 3D Trans Shade 7
MiYO 3D Trans Shade 8
MiYO 3D Trans Shade 9
MiYO 3D Trans Shade 10
MiYO 3D Trans Shade 11

VITA® classic refers to the VITA classical A1-D4® shade guide. VITA®, VITA classical A1-D4®, 3D-MASTER® and the related shade designations, are registered trademarks of Vita Zahnfabrik H. Rauter GmbH & Co. KG. MiYO®3D is not affiliated with or endorsed by VITA®.